

MEDICLEAN disinfectant for hospital use



Presentation:

100 tablets in pots of white polyethylene.

Use recommendations:

It is a product developed for use in hospitals, laboratories, dental centers, emergency rooms, dialysis centers and blood transfusions and other medical institutions in general.

Applications:

The Mediclean tablets were developed for use in disinfecting surfaces at risk of infection, quickly, accurately and safely.

Shelflife:

3 years

Advantages of this presentation:

- ◆ Easily metered, stable, accurate and safe concentrations tablets.
- ◆ Disinfection procedure operates under validating ppm chlorine (mg / liter).
- ◆ Allows use in conjunction with anionic detergents, soaps and degreasers for cleaning and disinfecting liquids simultaneously.

TABLETS OF 2,5 g. NaDCC (each tablet 1550 ppm free chlorine per liter of water)



LEVELS	Nº TABLETS	LITERS QUANTITY	ppm AVAILABLE CHLORINE
High risk of infection	6	1	9.500
	3	1	4.650
Medium risk of infection	1	1	1550
	1	2	775
Low risk of infection	1	5	310
	1	8	200



Ezpeleta 1277, office 1 | (B1640FXI) | Martínez | Buenos Aires | Argentina
Ph./Fax: +54 11 4798-0708 | e-mail: info@laboratoriopyam.com | web: laboratoriopyam.com
f /laboratoriopyam